irit®

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**IR-3602**

**Main functions**

- Takes off the effect of weariness.
- Stimulates blood circulation.
- Relaxes and relieves muscle stress.

**Safety Instructions**

- **Attention: this appliance is not intended for professional use.**
- Read the manual carefully before the first use.
- DO NOT immerse the appliance in water or other liquids, It may cause serious damage.
- Keep the massager out of the reach of children.
- Never leave the massager unattended, when it is switched on.
- Never try to repair the appliance by yourself.
- Keep the appliance away from the high temperature sources, direct sunlight and moisture.
- Avoid big pressure on the appliance, it may cause some damage.
- **Warning:** If the massager using caused any cramps, consult with your doctor before the next use.
- It is also recommended to stop using the massager, if you have any skin damage or dermatitis.
- The appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities or lack of experience and knowledge, unless they have been given supervision or instruction concerning the use of the appliance by a person responsible for their safety.

**Usage**

- Insert 3xAAA batteries (not included) into the batteries compartment.
- Either connect the massager with USB cable.
- Press ON button located on top of unit to start the massager.

**Cleaning**

- Use soft fabric to clean the appliance.

The importer retains his right to make changes in appearance and completeness of the item without preliminary notification to the customer.

For household usage only.

Do not leave unattended!

The date of manufacture is marked on the packaging

Service life is 2 years

Made in China

**SPECIFICATION:**

Power: 2 W

Powered by 3xAAA batteries (not included), usb cable

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